



**IHE Integration Statement for
I-SoftView 6.0**

September 25, 2004

I-SoftView is sold by Toshiba Medical Systems Europe
under the marketing name ReportDirect



Cedara I-SoftView IHE Integration Statement		Date
		<i>September 25, 2004</i>
Vendor	Product Name	Version
<i>Cedara Software Corp.</i>	<i>Cedara I-SoftView (Report)</i>	<i>6.0</i>
<p><i>This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options specified below:</i></p>		
Integration Profiles Implemented	Actors Implemented	Options Implemented
<ul style="list-style-type: none"> • Schedule Workflow (SWF) • Consistent Presentation of Images (CPI) • Access to Radiology Information (ARI) • Patient Information Reconciliation (PIR) • Key Image Note (KIN) • Evidence Documents (ED) 	Image Display	None
<p><i>Links to Standards Conformance Statements of the Implementation</i></p>		
HL7		
DICOM		
<p><i>Link to further information on this product:</i></p>	<p>http://www.cedara.com/products/diagnostic_and_review_workstations.htm</p>	